# Hotplates & Stirrers

Instructions for use

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	<b>Model No.</b>
Hotplates	SB300
	SB500
	CB300
	CB500
Digital Hotplates	SD300
	SD500
Magnetic Stirrer	SB301
Stirrer Hotplates	SB302
	CB302
Infrared Hotplate	CR300
Infrared Hotplate Stirrer	CR302

## Before Use

**Thank you for purchasing this Bibby Scientific product. To get the best performance from the equipment, and for your own safety, please read these instructions carefully before use.**

**If the equipment is not used in the manner described in this manual the protection provided may be impaired.**

This equipment is designed to operate under the following conditions:-

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuation not exceeding 10%
- ❖ Overvoltage category II IEC60364-4-443
- ❖ Pollution degree 2 IEC664
- ❖ Use with a minimum distance all around of 200mm from walls or other items

## Electrical Installation

Before connection please read and understand this instruction manual and ensure that the line supply corresponds to that shown on the rating plate. All models require a supply rated at: 230V, 50Hz, ~, single phase.

Power consumption is:

Model	Power	Model	Power
SB300	600W	CB300	1200W
SB301	50W	CB302	1250W
SB302	650W	CB500	2250W
SB500	1500W	CR300	900W
SD300	600W	CR302	950W
SD500	1500W		

**THIS INSTRUMENT MUST BE EARTHED**

All models are supplied with two mains leads fitted with IEC plugs for connection to the instrument. One has a U.K. 3 pin plug and the other has a 2 pin "Shuko" plug for connection to the mains. Choose the lead appropriate for your electrical installation and discard the other.

Should neither lead be suitable, take the lead with the U.K. plug and replace the plug with a suitable alternative. This involves cutting off the moulded plug, preparing the cable and connecting to the rewirable plug in accordance with its instructions.

**IT IS IMPORTANT THAT THIS OPERATION SHOULD ONLY BE UNDERTAKEN BY A QUALIFIED ELECTRICIAN**

NOTE: Refer to the equipment's rating plate to ensure that the plug and fusing are suitable for the voltage and wattage stated

The wires in the mains cable are coloured as follows:

BROWN	- LIVE
BLUE	- NEUTRAL
GREEN/YELLOW	- EARTH

The instruments are fitted with an IEC socket at the rear of the instrument for connection of the mains lead. The appropriate mains lead should be connected BEFORE connection to the mains supply.

Should the mains lead need replacement a cable of 1mm<sup>2</sup> of harmonised code H05W-F connected to an IEC320 plug should be used.

N.B. The UK mains lead is protected by a 10A fuse mounted in the plug top.

### **IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN**

When connected to the mains electricity supply a green light on the front panel illuminates, except for the digital SD300 and SD500 models where the display will illuminate showing the word "OFF".

## **Safety Advice Before Use**

- ❖ If the mains supply is interrupted the units will restart with its current settings on the restoration of the electricity supply.
- ❖ Do not use an instrument with a ceramic top if the top plate is scratched, chipped or shows signs of chemical etching.
- ❖ Position the instrument on a firm level surface away from any heat sensitive or flammable material.
- ❖ Do not use the units with flammable liquids or in a hazardous atmosphere.
- ❖ There is a danger of liquid spills if vessels are stirred too vigorously. Always build stirrer speed slowly and never stir faster than necessary.
- ❖ Never move or carry the instrument until it has been switched off and allowed to cool.
- ❖ Never move or carry the unit with containers on the top plate or while the instrument is connected to the mains electricity supply.

## **Hot Surfaces**

**WARNING: The top surface of the instrument may be HOT. Do not leave heaters switched on when not in use.**

For user safety all heated models are fitted with a safety light to advise when the top is hot. This light will flash red on the front panel when the plate temperature exceeds 50°C and will continue to flash even if the heater is switched off so long as the instrument remains connected to the electricity supply.

## **Heating**

**Models SB300, SB302, SB500, CB300, CB302, CB500:**

The knob labelled "Heat" on the front panel controls the heat output and hence the plate temperature.

The knob is graduated from 1 - 9. Increasing the value increases the heat output.

Turn the hotplate control knob clockwise to the desired setting. The amber load indicator light will illuminate and it will begin to cycle on and off when the desired temperature is reached and the hotplate begins to control.

**Models SD300 and SD500:**

When connected to the mains electricity supply the digital display will illuminate and show the word OFF.

Switch the unit on by pressing the left hand one of the 2 buttons located to the left of the digital display. The display will now show the actual plate temperature in degrees Celsius.

Press and hold the "press to set" button, the right hand one of the 2 buttons located to the left of the digital display. The display will now show the set temperature.

The set temperature may be adjusted by turning the knob located to the right of the digital display whilst continuing to hold in the "press to set" button. When the desired setting is displayed release the "press to set" button. The display will revert to show the actual plate

temperature and the instrument will begin to heat.

The set temperature can be observed at any time by pressing the "press to set" button.

Turning the adjusting knob without pressing the "press to set" switch will have no effect on the hotplate setting.

#### **Models CR300 and CR302:**

The knob labelled "Heat" on the front panel controls the heat output and hence the plate temperature. When the heat is switched on the amber pilot light illuminates.

The knob is graduated from 1 - 9. Increasing the value increases the heat output.

In operation only the area defined by the pattern will become fully heated. The outer edges remain cooler for safety considerations.

Vessels larger than the described area should not be used. Failure to observe this could give rise to overheating of the outer parts of the instrument.

**NOTE: To prevent overheating of the infrared heater a temperature limiter is fitted which switches off the heater should the temperature exceed 580°C.**

## **Stirring**

#### **Models SB301, SB302, CB302 and CR302:**

The knob labelled "Stir" on the front panel controls the stirrer speed.

It is labelled with an arbitrary numeric scale, 1 - 9 and turning the knob to a higher setting increases stirrer speed.

When the stirrer is switched on the amber pilot light illuminates.

## **Cleaning and Care**

Before attempting cleaning disconnect from the mains electricity supply and ensure that the top is cool.

Periodically clean the body of the instrument using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents.

**Ceramic top models:** The ceramic top is very easy to keep clean. A damp cloth will normally remove most types of contamination. For more difficult stains a domestic cream cleanser is recommended. Cleaning is made easier if spillages are attended to promptly. In any case, spillages of alkali, phosphoric acid and hydrofluoric acid **MUST** be removed immediately as these chemicals can attack and damage the glass ceramic. Ensure that the appropriate safety precautions are observed.

During cleaning and general operation take care not to scratch the surface of the top plate as this could result in subsequent thermal breakage.

**WARNING: A ceramic top which is scratched, chipped, chemically etched or damaged must not be used.**

**Metal top models:** The metal top plates should be cleaned using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents. Cleaning is made easier if spillages are attended to promptly. In any case, spillages of acids and alkalis **MUST** be removed immediately as these chemicals can attack and damage the aluminium alloy. Ensure that the appropriate safety precautions are observed.

## Servicing and Repair

This product range does not require any routine servicing.

**Note: There are no internal user replaceable parts**

In the event of product failure any repair or replacement of parts should only be undertaken by suitably qualified personnel. For advice or to receive a service manual please contact the Technical Service Department of Bibby Scientific Ltd. quoting model and serial number.

Only spare parts supplied or specified by Bibby Scientific Ltd. or its agents should be used. Fitting of non-approved parts may affect the performance of the safety features designed into the instrument.

**Note: The magnetic stirrer drive utilises strong magnets.**

If in any doubt, please contact the Technical Service Department of Bibby Scientific Ltd. or the point of sale.

Bibby Scientific Ltd.  
Stone, Staffordshire ST15 0SA  
United Kingdom  
Tel: +44 (0) 1785 812121  
Fax: +44 (0) 1785 813748  
e-mail [info@bibby-scientific.com](mailto:info@bibby-scientific.com)  
[www.bibby-scientific.com](http://www.bibby-scientific.com)

## Warranty

Bibby Scientific Ltd warrants these instruments to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of three (3) years. In the event of a justified claim Bibby Scientific will replace any defective component or replace the unit free of charge.

This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by installation, adaptation, modification, fitting of non approved parts or repair by unauthorised personnel.

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These products meet the relevant EC harmonised standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in their vicinity will meet these standards

and we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, please contact the Technical Service Department of Bibby Scientific Ltd.

## Declaration of Conformity

### Large Hotplates

#### Models

**SB300 SB301 SB302 SB500  
CB300 CB302 CB500  
SD300 SD500 HT3 HT4**

These products comply with the requirements of the EU Directives listed below:

**89/336/EEC** Electromagnetic Compatibility Directive amended by 93/68/EEC.

**73/23/EEC** Low Voltage Directive amended by 93/68/EEC

Compliance with the requirements of these Directives is claimed by meeting the following standards:

EN 61326: 1997 + Amendments A1: 1998, EN 61000-4-2: 1995, EN 61000-4-4: 1995, EN 61000-4-3: 1996, EN 61000-4-5: 1995, EN 61000-4-6: 1996, EN 61000-4-11: 1994, EN 61000-3-2: 2001, EN 61000-3-3: 1995, EN 61010-1: 2001, EN 61010-2-10: 1994, EN 61010-2-51: 1995.

#### Compliance Certificates and Full Reports.

Ref: RETS0505/A/1, RETS0505/A/2, RETS0505/A/3, RETS0505/A/4, RETS0505/A/5, RETS0505/A/6, RETS0505/A/7, RETS0505/A/8,

From "Epsilon Technical Services" an Independent Accredited Test House Showing Compliance to the above Standards are available on request.

#### **CE Mark Affixed 02**

Signed:

( Mr D E Hicks )

Date:

25/7/05

Authority:

Technical / Development Manager

For

Barloworld Scientific Ltd - Stone - Staffs - ST15 0SA - UK  
Tel: +44 (0) 1785 812121 - Fax +44 (0) 1785 813748

# Declaration of Conformity

## Infra-red Hotplate & Infra-red Hotplate Stirrer

### Models

CR300, CR302

These products comply with the requirements of the EU Directives listed below:

**89/336/EEC**                      Electromagnetic Compatibility Directive amended by 93/68/EEC.

**73/23/EEC**                      Low Voltage Directive amended by 93/68/EEC.

Compliance with the requirements of these Directives is claimed by meeting the following standards:

EN 61326: 1997 + Amendments A1 & A2.  
(Electrical Equipment for Measurement, Control and Laboratory use).

EN 61010-1: 2001.  
(Safety Requirements Electrical Equipment for Measurement, Control and Laboratory use).

EN 61010-2-010: 2003.  
(Particular requirements for laboratory equipment for heating of materials).

### Compliance Certificates and Full Reports

Ref: RETS(E)1002/A/1, RETS(L)1002/A/2, from 'Epsilon Technical Services'  
(an Independent Accredited Test House) showing compliance to the above standards, are available on request.

**CE mark affixed '04.**

Signed:  (Mr D E Hicks)

Date: 26/7/05

Authority:      Technical / Development Manager

For

Barloworld Scientific Ltd - Stone - Staffs - ST15 0SA - UK  
Tel: +44 (0) 1785 812121 - Fax +44 (0) 1785 813748

# INSPECTION REPORT

MODEL \_\_\_\_\_

SERIAL No. \_\_\_\_\_


## ELECTRICAL SAFETY

- 1. Earth continuity
- 2. Insulation
- 3. Flash test

## FUNCTIONAL

- 1. Indicators
- 2. Heater/stirrer control
- 3. Visual acceptance

QUALITY CONTROL INSPECTOR

### **Bibby Scientific France SAS**

ZI du Rocher Vert - BP 79  
77793 Nemours Cedex  
France  
Tel: +33 1 64 45 13 13  
Fax: +33 1 64 45 13 00  
e-mail: bsf@bibby-scientific.fr  
www.bibby-scientific.com

### **Bibby Scientific Italia Srl**

Via Alcide de Gasperi 56  
20077 Riozzo di Cerro al Lambro  
Milano Italia  
Tel: +39 (0)2 98230679  
Fax: +39 (0)2 98230211  
e-mail: marketing@bibby-scientific.it  
www.bibby-scientific.it

### **Bibby Scientific Ltd**

Beacon Road Stone  
Staffordshire ST15 0SA  
United Kingdom  
Tel: +44 (0)1785 812121  
Fax: +44 (0)1785 813748  
email: info@bibby-scientific.com  
www.bibby-scientific.com

### **Bibby Scientific US Ltd**

3 Terri Lane Suite 10  
Burlington NJ 08016  
USA  
Tel: 800-225-9243  
Fax: 609-589-2571  
www.bibby-scientific.com

### **Afora S.A.**

Calle Aribau 240  
08006 Barcelona Spain  
Tel: +343 93-306 98 00  
Fax: +343 93-306 98 23  
e-mail: marketing@afora.com  
www.afora.com