# Overhead stirrer

SS10 SS20 SS30


Instructions for use

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
**Figure 1**

**Thank you for purchasing this Bibby Scientific product. To get the best performance from the equipment, and for your personal safety, please read these instructions carefully before use.**

This equipment is designed to operate under the following conditions:-

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuation not exceeding 10%
- ❖ Overvoltage category II IEC60364-4-443
- ❖ Pollution degree 2 IEC664

If the equipment is not used in the manner described in this manual the protection provided by the equipment may be impaired.

 **WARNING!** This Overhead Stirrer is not suitable for use in inflammable or explosive atmospheres. It should not be used to stir any inflammable liquids unless they are in a suitable sealed vessel where there is no possibility of gases or vapours coming into contact with the Overhead Stirrer body.

## Electrical Installation

 **THIS EQUIPMENT MUST BE EARTHED**

**Before connection please ensure that the mains power supply corresponds to that shown on the rating plate.**

All models are supplied with a mains lead fitted with an IEC plug to fit the IEC socket on the short flying lead emerging from the Stirrer mounting bar. This lead should be connected to the flying lead BEFORE connection to the mains supply. Should the mains lead require replacement a cable of 1mm<sup>2</sup> of harmonised code H05W-F connected to an IEC 320 plug should be used.

Most Stuart equipment is supplied with the mains lead fitted with a moulded plug for connection to the electricity supply. If this is not the case it should be fitted with a plug approved for the country in which it is to be used, noting that the mains lead has cables coloured according to the following code:

Blue-	Neutral
Brown-	Live
Green/yellow-	Earth

If required the plug should be fitted with a 10A fuse to protect the cable. All models are protected by two fuses located in the back panel. These conform to IEC127 and are rated as follows:

Model	Power	Fuse
SS10	50W	T1A
SS20 & SS30	80W	T1A

**IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN**

## General Description

The Stuart Overhead Stirrers are designed for effective mixing of large volumes or viscous liquids. All models have an easy to use keyless chuck, an enclosed mains cable designed to be kept away from moving parts, and overload protection for maximum safety.

**SS10 General Purpose Stirrer:** is designed for use in aqueous and low viscosity liquids to a maximum volume of 15 litres.

**SS20 High Performance Stirrer:** is designed for use in low to medium viscosity liquids to a maximum volume of 25 litres. It has a hollow spindle to enable long stirrer shafts to be used or for raising the paddle out of vessels easily.

**SS30 Dual Torque Stirrer:** is designed for use in low to high viscosity liquids to a maximum volume of 40 litres. It has the same hollow spindle as the SS20, and two modes of operation - mode 1 gives high torque at lower speeds, mode 2 higher speed at lower torque.

Viscosity	mPa-s	Example
Very low	0 - 100	Water to thin oil
Low	100 - 1000	Thin to thick oil
Medium	1000 - 10,000	Thick oil
High	10,000 - 50,000	Thick oils to honey

## Safety Advice Before Use

- ❖ Never operate the Stirrer without first clamping the mounting rod securely. It is NOT intended to be hand held. In use the stability of the Overhead Stirrer is dependant on the security of the mounting assembly. The operator should ensure that the mounting is adequate for the intended application. (See the Accessories Section for stand and stirrer paddles available).
- ❖ Clamp the stirrer shaft securely in the chuck and lock before operating the Stirrer.
- ❖ Do not operate models SS20 or SS30 with a stirrer shaft protruding through the top of the Overhead Stirrer without first removing the top grommet.
- ❖ Do not mix flammable liquids or use the equipment in hazardous atmospheres.
- ❖ In case of mains interrupting, a fault or mechanical failure, the unit will continue to operate on removal of fault.

- ❖ Mechanical energy can lead to breakage of glass vessels. Use with care.

**WARNING!** The Overhead Stirrer will automatically restart after a power failure. It may also restart in the event of a fault or mechanical interruption (stall).

## Operation

### Fixing a Stirrer Paddle

The chuck is designed to grip a stirrer shaft up to a maximum of 13mm diameter. Loosen the jaws of the chuck by turning clockwise. Place a stirrer shaft in the chuck and tighten by turning anti clockwise. The chuck only requires tightening by hand.

The chuck protective cover is designed for use with stirrer shafts up to 8mm diameter only. It grips the top of the shaft preventing ingress of corrosive materials. If a stirrer shaft with a diameter greater than 8mm is used, remove the protective cover before fixing the shaft in the chuck.

### Hollow Spindle (SS20 & SS30 only)

Models SS20 & SS30 have a hollow spindle enabling a stirrer shaft of diameter up to 8mm to be passed through the Stirrer body and adjusted from the top. Open the chuck and remove the top grommet. Feed the shaft through the chuck and into the Stirrer until it protrudes through the top. Tighten the chuck as described above.

### Using the Stirrer (SS10, SS20 & SS30)

Clamp the Overhead Stirrer securely to a stand and fix a stirrer paddle in the chuck (as described above). Check that the Stirrer is switched OFF before inserting the mains cable into the short flying lead and connecting to the mains electricity supply. When the supply is available from the mains, a green light will illuminate on the front of the Stirrer.

### Selecting the Stirring Speed (SS10 & SS20)

Switch the Stirrer ON by the mains on / off switch. A yellow light will illuminate on the front of the Stirrer. Rotate the speed control knob to select the required speed of rotation. A soft start facility is built into the Stirrer so there will be a slight delay between selection and rotation.

### Selecting the Speed / Torque mode (SS30 only)

Switch the SS30 Stirrer on by the mains on / off switch. A yellow light will illuminate on the front of the Stirrer. Rotate the speed control knob to select the required mode:

- Clockwise selects the high speed lower torque mode (max. speed/ torque 2000rpm / 27Ncm)
- Anticlockwise selects the low speed higher torque mode (max. speed/torque 500rpm / 90Ncm).

A soft start facility is built into the Stirrer so there will be a slight delay between selection and rotation and when switching between torque modes.

### Overload Protection (SS10, SS20 & SS30)

If the load on the Overhead Stirrer is too great a red warning light on the front of the Stirrer will illuminate and the power to the motor will be limited to protect the Stirrer from damage. The Stirrer will continue to run at the limited speed until the load is reduced.

If there is a danger of overheating the power to the motor will be further reduced. If the danger of overheating becomes severe the Stirrer will stop. The unit must be switched OFF by the mains on / off switch and allowed to cool before it can be used again.

## Maintenance & Servicing

**WARNING:** Ensure the unit is disconnected from the mains electricity supply before attempting maintenance or servicing.

Periodically clean the instrument using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents.

**Any repairs or replacement of parts MUST be undertaken by suitably qualified personnel.**

## Spares

The following spares and accessories are available from your laboratory dealer:

Description	Catalogue Number
Fuse, T1A	ESS20C1(S)
Top case grommet	USS20C5
Chuck protective cover	USS2056

### For complete paddle, order rod plus head(s)

Paddle rod, 350mm	SS10/5
Paddle rod, 550mm	SS10/6
Small paddle head	SS10/10
Large paddle head	SS10/11
Large propeller head (2 blade)	SS10/12
Small propeller head (4 blade)	SS10/13
Anchor paddle head	SS10/14
Stand, General Purpose (inc. base, rod & bosshead)	SS10/1
Stand, Heavy Duty (inc. base, rod & bosshead)	SS10/2

For a comprehensive list of parts and a Service Manual required by service engineers conducting internal repairs, please contact the Technical Service Department of Bibby Scientific Ltd. quoting both the model and serial number.

Only spare parts supplied or specified by Bibby Scientific Ltd. or its agents should be used. Fitting of non-approved parts may affect the performance and safety features designed into the instrument.

## Technical Data

	SS10	SS20	SS30	
			Low torque	High torque
Speed range	100 - 2000rpm	100 - 2000rpm	100 - 2000rpm	50 - 500rpm
Litre capacity	15L	25L	25L	40L
Max. viscosity	10,000mPas	20,000mPas	20,000mPas	50,000mPas
Max. torque	15Ncm	27Ncm	27Ncm	90Ncm
Input power	50W	80W	80W	
Chuck range	1.5 to 13mm	1.5 to 13mm	1.5 to 13mm	
Hollow spindle	-	8.2mm dia.	8.2mm dia.	
Motor type	DC permanent magnet	DC permanent magnet	DC permanent magnet	

### Dimensions

- Stirrer (wxdxh)	85 x 175 x 230mm	85 x 195 x 230mm	85 x 195 x 230mm
- Mounting bar	12.7 dia. x 240mm	12.7 dia. x 240mm	12.7 dia. x 240mm
Net weight	1.9Kg	2.1Kg	2.5Kg

If in any doubt, please contact the Technical Dept. of Bibby Scientific or the point of sale.

Bibby Scientific Ltd.

Beacon Road Stone  
Staffordshire ST15 0SA  
United Kingdom

Tel: +44 (0)1785 812121

Fax: +44 (0)1785 813748

e-mail: info@bibby-scientific.com

www.bibby-scientific.com

## Warranty

Bibby Scientific Ltd warrants this instrument to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of three (3) years. In the event of a justified claim, Bibby Scientific will replace any defective component or replace the unit free of charge.

This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by installation, adaptation, modification, fitting of non-approved parts or repair by unauthorized personnel.


Figure 1

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# Notes

These products meet the relevant EC harmonised standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in their vicinity will meet these standards and we cannot guarantee that interference will not occur

in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, please contact the Technical Service Department of Bibby Scientific Ltd.

## ***Declaration of Conformity***

---

### **Overhead Stirrers**

#### **Model SS10, SS20 & SS30**

These products comply with the requirements of the EU Directives listed below:

**89/336/EEC** Electromagnetic Compatibility Directive, amended by 92/31/EEC and 93/68/EEC.

**73/23/EEC** Low Voltage Directive amended by 93/68/EEC

Compliance with the requirements of these Directives is claimed by meeting the following standards:

EN 61326: 1997 + Amendments A1: 1998, EN 61000-4-2: 1995, EN 61000-4-4: 1995, EN 61000-4-5: 1995, EN 61000-4-6: 1996, EN 61000-4-9: 1993, EN 61000-4-11: 1994, EN 61000-3-2: 1995, EN 61000-3-3: 1995, EN 61010-1: 1993 + Amendments A1: 1992 & A2: 1995, EN 61010-2-51: 1995.

#### Compliance Certificates and Full Reports.

Ref: RETS0287/47 and RETS0287/48

From "Epsilon Technical Services", an Independent Accredited Test House, showing compliance to the above standards, are available on request.

**CE mark affixed '01.**

Signed:  (Mr D. E. Hicks)

Date: 18/10/05

Authority: Technical/Development Manager

for

Barloworld Scientific Ltd - Stone - Staffs - ST15 0SA - UK  
Tel: +44 (0) 1785 812121 - Fax +44 (0) 1785 813748

# INSPECTION REPORT

**MODEL** \_\_\_\_\_


## ELECTRICAL SAFETY

- 1. Earth continuity 
- 2. Insulation 
- 3. Flash test 

## FUNCTIONAL

- 1. Indicators 
- 2. Speed control 
- 3. Visual acceptance 

**QUALITY CONTROL INSPECTOR**

### **Bibby Scientific France SAS**

ZI du Rocher Vert - BP 79  
77793 Nemours Cedex  
France  
Tel: +33 1 64 45 13 13  
Fax: +33 1 64 45 13 00  
e-mail: [bsf@bibby-scientific.fr](mailto:bsf@bibby-scientific.fr)  
[www.bibby-scientific.com](http://www.bibby-scientific.com)

### **Bibby Scientific Italia Srl**

Via Alcide de Gasperi 56  
20077 Riozzo di Cerro al Lambro  
Milano Italia  
Tel: +39 (0)2 98230679  
Fax: +39 (0)2 98230211  
e-mail: [marketing@bibby-scientific.it](mailto:marketing@bibby-scientific.it)  
[www.bibby-scientific.it](http://www.bibby-scientific.it)

### **Bibby Scientific Ltd**

Beacon Road Stone  
Staffordshire ST15 0SA  
United Kingdom  
Tel: +44 (0)1785 812121  
Fax: +44 (0)1785 813748  
e-mail: [info@bibby-scientific.com](mailto:info@bibby-scientific.com)  
[www.bibby-scientific.com](http://www.bibby-scientific.com)

### **Bibby Scientific US Ltd**

3 Terri Lane Suite 10  
Burlington NJ 08016  
USA  
Tel: 800-225-9243  
Fax: 609-589-2571  
[www.bibby-scientific.com](http://www.bibby-scientific.com)

### **Bibby Scientific (Asia) Ltd**

Room 607 Yen Sheng Centre  
64 Hoi Yuen Road Kwun Tong  
Kowloon Hong Kong  
Tel: +852 3583 1581  
Fax: +852 3583 1580  
e-mail: [bibby@bibby-scientificasia.com](mailto:bibby@bibby-scientificasia.com)  
[www.bibby-scientificasia.com](http://www.bibby-scientificasia.com)